

HOW TO

Download the Full Viewer

- 1** Click the Account link on the main toolbar.
(gear icon near the top right).



- 2** Click the Image Viewer tab.

Account Settings


Image Viewer

- 3** Click the Download Full Viewer button.


DOWNLOAD FULL VIEWER



- A. Save the file when prompted.



- B. Display the downloads by clicking on the down arrow at the far top right of the screen.




- C. Click on eradsetup.exe

OR




- A. Click to "Run" the file when prompted to run setup wizard.




OR



- A. Right click "eradsetup.exe" in the lower left corner of the screen.
B. Click "open" to run setup wizard.



- C. Follow the setup wizard, accepting the License Agreement and accepting all default options. Note: If the user of this PC does not have admin privileges, you will receive an error at the start of the wizard. The office's IT will need to be contacted. The WebTeam can assist with this.



- D. A reboot may be required if the PC is running Windows 7 and if an error is received on launch.
E. Run the PACS Viewer for the first time so all necessary updates can occur.